



1. Product name

Southpaw Acrobat Swing

2. Product code

21009065

3. Colour

As shown



4. Short description

The Acrobat Swing is a unique and versatile tool specifically designed for sensory stimulation, motor training, and relaxation. Thanks to its different suspension options, the swing can be used in various ways, tailored to the needs of the user.

5. Contents

- 1x Southpaw Acrobat Swing
 - 1x maintenance checklist
 - 1x safety snap instruction sheet
- WORKING LOAD: 100 lbs. per layer / 400 lbs. total.

6. How to use

The Southpaw Acrobat Swing can be suspended in several ways:

1. All four corners can be connected to a single suspension point, like a net swing. *Note: due to the stretchiness of the fabric, a high ceiling is required for this configuration.*
2. It can be suspended like a hammock from two points.
3. It can be attached with two points on one end and one on the other, in a triangle configuration.
4. Each corner can be attached to its own eyebolt.

Important: Do not use the Acrobat Swing while carrying objects that could damage the fabric, such as a pencil or a comb. Also, do not sit directly on the stitched hem, as this may cause the seams to tear.

Warning:

Failure to properly install a forged steel eye bolt (1/2" or 5/8") may result in serious injury or death. If you have any questions regarding the installation of the forged steel eye bolt, please contact us.

You can use the **Southpaw Acrobat Swing** in a variety of ways, for different therapeutic goals. To suspend the swing in one of the configurations mentioned above, simply attach a **Southpaw Suspension System – Safety Snap (13799652)** to each webbing loop at the corners of the swing, and then connect the snap to your suspension point.

If you want to attach more than one corner to a single eyebolt, connect those webbing loops to the same safety snap. A **Southpaw Suspension System – Height Adjuster (13626652)** may be required to adjust one or more corners to the proper height for the activity or client.

When all four corners of the swing are attached to a single suspension point, a **Southpaw Suspension System – Rotational Device (13738652)** must always be used.





Note: The user may need assistance when getting in and out of the swing.

Space Requirements:

Always maintain a minimum of **1.8 meters (6 feet)** of clear space between the swing and walls or other obstacles.

With proper supervision, however, the swing can also be positioned in combination with other equipment to allow motor planning activities where the client moves from one apparatus into the swing, or from the Acrobat Swing onto another swing.

The choice of activity should always be based on the experience of the therapist or parent, and their ability to provide a challenging yet safe activity.

7. Maintenance

Monthly checks (record on checklist):

All the following wear points should be checked on a monthly basis and recorded on the checklist:

- Check the webbing for wear, especially at the point where it rubs against the Southpaw Suspension System – Safety Snap. If wear is found: discontinue use until repaired.
- Check the corners of the Acrobat Swing. These areas are critical. If the stitching is coming loose, or the Lycra is torn or pulling out of the corner: discontinue use until repaired.
- Check the seams for loose or missing stitches. If damage is found: discontinue use until repaired.
- Refer to the other instruction sheets for the proper maintenance of related items, such as the Southpaw Suspension System – Safety Snap and the Southpaw Suspension System – Rotational Device.

Matting – Precautions

- Always use appropriate crash mats.
- Place crash mats and wall mats with all suspended equipment.

Storage

- Store away from direct sunlight
- Keep away from chemicals and cleaning agents
- Not accessible to clients without supervision
- Do not store on the edge (may damage the rim)

Cleaning

- Check with your infection control department before using these procedures.
- Check first for client allergies.
- Machine wash in warm water with mild detergent. Tumble dry on low heat.





8. Safety

The safety of your clients depends on carefully selecting activities that match — and challenge — their physical abilities.

Always use proper spotting and safety techniques. In addition, the layout of the room, correct installation, and regular maintenance of your equipment are critical.

Southpaw equipment is not intended for recreational or outdoor use. Under no circumstances should this equipment be used without supervision.

This document is for example purposes only and in no way replaces local, regional, national, or international safety standards. Please note: it is your responsibility to obtain all relevant safety information that applies to your facility.

For more information, you can visit our website www.nenko.com

